
BoneFlo®

Chirurgisches Saugerset / Surgical Suction Set

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Please read the Instructions for Use (IFU) carefully before using BoneFlo®

EN

Surgical suction set BoneFlo®

Item number: B2017M0

English (EN) ISO-639-1

Instructions for use (IFU)

REF


A. Product description

BoneFlo® is a commercial, bio-compatible surgical suction set made of polymethyl methacrylate (PMMA). This single-use product is suitable for the removal of small, detached tissue components (bone and soft-tissue fragments, bone marrow, blood, bodily fluids, clotted blood, etc.) or foreign bodies and was developed primarily for orthopaedic and trauma surgery. The suction set is connected to a receptacle via a tube before negative pressure is applied.

The suction tip (1) (a cannula angled at 20°) is connected to the system's inlet opening (2), while conventional commercial tube systems Ch25 can be connected to the conical outlet opening (3). BoneFlo® also contains a replaceable internal plastic filter (4) which traps hard tissue components, clotted blood and foreign bodies during the suction procedure. The special BoneFlo® design includes a screw cap (5) allowing replacement of the plastic filter. It is necessary to change the filter when it becomes clogged with tissue components or foreign bodies and suction performance decreases. Changing the filter involves opening the screw top in the outlet area and manually removing the filter. A manual tilting movement usually suffices here. New plastic filters are then applied, with the arrow marking visible on the filter pointing in the direction of the suction tip. The cap can then be screwed on with no further handling required thanks to the locking mechanism in the filter and the suction set is once again available for use on the surgical site. Cleaning the filter is not recommended.

The set includes a pre-assembled suction handle with a filter inset, an additional filter (loose) and two suction tips. Before use, the installed suction handle needs to be connected to the suction tip (inlet opening: 2) and to the surgical tube (outlet opening: 3). The stability of the plug-in connections needs to be controlled prior to application.

Product data sheet

Material specifications

Contains latex:	No
Suction set, filter, tip	PMMA, pyrogen free

Certification

CE marking:	Yes (CE0044)
Medical device classification:	Ila (EU)
Quality management certified in accordance with	ISO 13485
UMDNS code:	13-846

Sterilisation

Sterile:	Yes
Sterilised by:	γ irradiation
Storage life/Service life:	see label

This is a sterile single-use product. Do not re-use! Not suitable for re-sterilisation. Re-use may result in function being compromised (material fracture, distortion).

Packaging

Individual packaging
 CLEERPEEL® Double peel packaging
 Dimensions [mm]: L 380 x W 160 x H 38; qty [pcs]: 1


Cardboard goods container

Internal dimensions [mm]: L 430 x W 255 x H 195; qty [pcs]: 15

Gross weight [kg]: 1.58 kg; volume [m³]: 0.02138; MOQ [pcs]: 1

Instructions for double bag packaging

Package includes a double sterile barrier system

- 1) The outer sterile barrier (1. layer) is opened by non-sterile OR staff
- 2) The inner sterile barrier is opened in an aseptic manner and is presented to the sterile OR staff
- 3) The inner sterile barrier (2. layer) is received by the sterile OR staff
- 4) The inner sterile barrier (2. layer) is opened by the sterile OR staff and the set can be removed

Storage and handling

- No not store above +40°C or below -30°C
- Keep dry, protect from moisture
- Protect from direct exposure to light
- Do not use if packaging is damaged
- Please note the expiry date

B. Indications

Surgical procedures on areas of the body where relevant quantities of tissue are released (for example orthopaedic or trauma surgery)

C. Contraindications

Application of BoneFlo® in visceral surgery. Surgical procedures in anatomical regions requiring special suction systems, such as intradural or intracerebral operations.

D. Adverse effects

- Air-flow noise during suction procedure

E. Warning

1. Do not use if the packaging has already been opened or is damaged.
2. This is a single-use product. Not suitable for multiple use or re-sterilisation, as this may lead to application-related product failure with injury to patients and/or users.
3. BoneFlo® consists of a container and a filter which retains and collects tissue components and foreign bodies from the surgical area. Please read carefully through the data sheet and instructions for use for further products, if combining them with BoneFlo®. TissueFlow GmbH assumes no liability or responsibility for the safety or performance of other products if they are used together with BoneFlo®.
4. Patient incompatibility with the product should be considered before use.
5. Infectious waste such as contaminated filters must be disposed of professionally and in accordance with the applicable laws and the ordinances of the respective hospital or institution.
6. Serious incidents must be reported to TissueFlow GmbH or the respective national sales representative and responsible health authorities.
7. The application of BoneFlo® is only permitted if combined with CE-certified tubes and suction systems. Please, note the IFU of the relevant manufacturers.
8. Please note that the application of too high vacuum might damage tissue.

F. Combination with other products

- The connection of BoneFlo® with an CE-marked vacuum system with reservoir is performed by an approved OR connection tube with a universal conical taper size Ch25.
- BoneFlo® can be connected to all approved surgical vacuum systems.
- The recommended vacuum is at -500mbar (corresponding 50 kPa or -375mmHg).
- The manual plug-in connection of the outlet opening (3) to the conical taper of the OR tube needs to be checked by axial traction and if necessary, adapted to prevent from disconnection during application.

G. Necessary qualifications for use

BoneFlo® may only be used by or under the supervision of a qualified doctor, dentist or equivalently qualified individual.

For more information, please contact the responsible local sales representative or TissueFlow GmbH.

H. Group of patients

BoneFlo® removes body fluids and tissue fragments in musculoskeletal patients to allow the surgeon a sufficient overview to the situs during surgery. The application of BoneFlo® in infants or toddlers (< 6 years of age) is not recommended.

Subject to errors, printing errors and technical amendments.

Contact information

TissueFlow GmbH
Am Wiesental 22
D-45133 Essen, Germany

The TissueFlow logo is a trademark of TissueFlow GmbH

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BoneFlo® is a trademark of TissueFlow GmbH.
Daha fazla bilgi için lütfen sorumlu yerel satış temsilcisi veya TissueFlow GmbH ile iletişime geçiniz.
**TissueFlow** GmbH
TissueFlow logosu, TissueFlow GmbH'nin ticari markasıdır